

HOWARD™

MEDICAL



Redefining medical technology.

Why Howard Medical

With over 35 years of world-class manufacturing experience, we design our products and services to ensure that we meet your specific needs and provide you with a point-of-care solution that is just that – a solution. Our extensive portfolio includes ergonomically-designed medication carts and nurse charting solutions that are fully integrated with tablets, notebooks, LCD monitors, and the accessories required for optimum job performance. And we offer FREE 100% US-based customer service and support for the life of your equipment. Our innovative products are real-world solutions for today's health care personnel.



Our offerings include:

- Fully customizable **medication boxes** that can be custom-built to your exact specifications.
- A variety of articulating **arms** that allow you to position equipment right where you need it.
- A host of **accessories** such as halogen lights, Sharps containers, bar code scanners, and much more.
- A vast selection of **technology equipment** from notebooks and tablets to desktops, printers, and faxes.

For more information about us, our products, or services, please call us at 877.856.6441 or visit our web site at www.Howard-Medical.com

Buy online at www.Howard-Medical.com or call us at 877.856.6441.

HOWARD'S ERGONOMIC POINT-OF-CARE CARTS



Version 2

Version 1

HI-Mobile

Compact and light, the HI-Mobile allows you to add shelves and bins to give you optimum storage capacity. Available in two versions.



Medium

Small

HI-Med

Supporting the largest bins available (5.5" to 8.5"), the HI-Med can be configured to your exact specifications. This spacious cart comes in small, medium, or large.




HI-Care

Our revolutionary HI-Care cart leads the industry with features such as tool-less battery replacement and conductive casters that help control electrostatic discharge.

EXPERIENCE THE HOWARD ADVANTAGE

At Howard Medical, we strive to bring to market a custom-built cart that is easy to use, custom-configurable, and one that is a necessity in the marketplace. That's why when we begin designing our carts, we don't just pull a design out of a hat. We get out and speak with nurses and hospital personnel to get their feedback on features they would like to see on carts that would make their jobs easier. We understand that each health care facility is unique, and we believe that we have a custom-configurable cart that will meet the most demanding specifications. Our line of carts allows nurses and hospital personnel to perform their duties without worry so that their main focus remains on what matters most...their patients.

SPECIAL FEATURES


- 2-point, swivel-mount monitors allow for comfort when wearing bifocals/contacts
- 270 degree handle enables easy steering and maneuverability
- Near-silent mechanical lift raises and lowers work surface making it possible to work comfortably while sitting or standing
- Conductive casters help minimize electrostatic discharge issues
- Customized colors are available: 



Swivel monitor available on the HI-Care and the HI-Prodigy



270° handle available on the HI-Care and the HI-Prodigy



Mechanical lift available on the HI-Care

HI-Prodigy



Our newest addition, HI-Prodigy, offers the same great features as our HI-Care cart but with a custom-configurable drawer system.



Small

HI-Pinnacle

Medium



Available in small, medium, and large with custom configurable drawer designs, the HI-Pinnacle can be fully customized with bar code scanners, halogen lights, and other necessary equipment without compromising work area.

TRAINING & INTEGRATION

We know that training and integration play an important role when choosing a medcart provider. Our knowledgeable Professional Services staff can provide you with both so that your cart is customized to meet your needs and you feel comfortable operating your new equipment. And if you have existing accessories you want added to your new carts, we can do that, too! Just ask.




Conductive casters available on HI-Care, HI-Prodigy, and HI-Pinnacle



Don't worry...you have options!

7,256 different options to be exact. That's how many different ways you can configure our HI-Med and HI-Pinnacle carts. Everything from large bulk storage drawers to small patient bins, large cup holders to imagers and computing devices. We customize the carts, not only to your hospital, but to your specific ward or unit. You tell us what your requirements are and we do the rest... CONSIDER IT DONE!






DID YOU KNOW WE offer
 MORE THAN
190,000
 QUALITY PRODUCTS?

desktops



notebooks



tablets



HOWARD™

TECHNOLOGY SOLUTIONS

Howard Technology Solutions is a leading manufacturer and reseller of computer hardware and related equipment. With a 100% US-based workforce and rich portfolio of over 190,000 products - we define technology roadmaps.

www.Howard.com

cables



mice



displays



projectors



faxes



bar code scanners

