

Locking Options

Electronic Narcotics Drawer Locking-Available for Advanced Electronics Only (note: all narcotic drawers must be two tier drawers) *Narcotics locking drawers replace cassettes in pairs of tiers, starting from the top working down.

Keyed Narcotics Drawer Locking

*Any two-tier drawer may have optional keyed narcotics lock installed.

Advanced Electronic Locking Features

- UP to 500 programmable user ID's
- PC Interface with Med Display Software
- LCD touch screen access
- Local and Remote Management of Hub settings and user ID's with Med Manager
- Automatic Relock with adjustable timeout
- Keyed override lock
- Drawer sense and audit capability
- Electronic or keyed narcotic drawer upgrades

Basic Electronic Locking Features

- 1 Clinician and 1 Pharmacy ID
- Keypad access
- Automatic relock with adjustable timeout
- Keyed override lock
- Keyed narcotics drawer upgrades

Push Button Combination Locking Features

- Keyless entry
- Keyed narcotics drawer upgrades

Keyed Locking

- Simple and secure
- Keyed narcotics drawer upgrades